

Eye Surgery Stretcher and Surgistool® Chair



Eye Surgery Stretcher

Access. Simple. Functionality.

Stryker's Eye Surgery Stretcher offers maximum headend clearance, giving surgeons ample leg room to maneuver while performing delicate eye procedures. The reliable work surface and mobile design combine to help with surgical throughput.



Enhanced head end clearance

Maximizes the surgeon's access for either a temporal or a superior approach to surgery.



Optional contoured or flathead piece comfortably supports patient during procedures.

Retractable fifth wheel helps improve mobility.

Metal tooth and groove braking system prevents swivel and roll.

Option of two- or three-position fold-down siderails.

Multiple bumper color options allow you to customize the stretcher for identification purposes.



Navy Gray Red

Surgistool Chair

The ideal surgeon's stool

Surgistool Chair offers customizable seat configurations with fully adjustable backrest and armrests provide optimal support to meet the surgeon's specific positioning preferences. Superior pedal access and a stable, yet mobile, design helps to facilitate clinical performance.



Three spindle base provides superior control and pedal access.



Optional padded armrests allows versatile positioning and stability.

Multiple seat configurations offer personalized comfort.



Teardrop seat



Wedge seat



Saddle seat
(optional)

Eye Surgery Stretcher and Surgistool Chair



Standard features: Eye Surgery Stretcher

- 3" (8 cm) Enhanced Comfort Mattress
- 8" (20 cm) Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake and steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

Standard features: Surgistool Chair

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Surgistool, Stryker. All other trademarks are trademarks of their respective owners or holder.

Copyright © 2017 Stryker
Mkt Lit-63 01262005 Rev H.2

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

patienthandling.stryker.com

Specifications

	Eye Surgery Stretcher	Surgistool Chair
Model number	1089	830
Overall length	90" (228.6 cm)	
Overall width	31.5" (80 cm)	
Weight capacity	500 lb. (225 kg)	300 lb. (136 kg)
Height range		
High	34" (86.4 cm) (to litter top)	28.5" (72 cm) / 30.5" (77 cm) (with riser kit)
Low	22.25" (56.5 cm) (to litter top)	21" (53 cm) / 23" (58 cm) (with riser kit)
Seat cushion style		Teardrop / wedge / saddle
Litter positioning		
Backrest	0°– 90°	
Knee gatch	0°– 30°	
Trend./Reverse Trend.	± 18°	
Patient surface	26" x 87" (66 x 221 cm)	
Siderails	13" x 55" (33 x 139.5 cm)	
Caster diameter	8" (20 cm)	

Stryker reserves the right to change specifications without notice.

Optional features: Eye Surgery Stretcher

- 4" (10 cm) Ultra Comfort Mattress
- Pre/post-op shoulder pads assembly (crank backrest only)
- Air delivery system/drape support
- Crank backrest
- Crank knee gatch
- Defibrillator tray/foot extender
- Temporal or superior wristrest
- Footboard/chartholder
- Three-sided hydraulic controls
- Three-sided brake and steer pedals
- Integrated patient transfer/arm board
- IV caddy
- Restraint straps
- Serving/instrument tray
- Siderail pads
- Surgery accessory rail
- Three-position siderails
- Upright oxygen bottle holder
- Oxygen bottle retaining ring
- Serving/instrument tray holder
- Head pads - Ultra Comfort
 - 3" (8 cm), flat
 - 3" (8 cm), concave
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable
- Extended warranty available

Optional features: Surgistool Chair

- 2" (5 cm) riser kit
- Saddle seat
- Multi-position armrest